



NIRS-SPM Instruction using Shimadzu Data




2009/12/07 Akihiro Ishikawa

Main Panel




Convert




Convert

We support to read our data format.




Shimadzu_FT_Right_4x4x2.TXT




For example...
Shimadzu_FT_Right_4x4x2.mat

Time Series Analysis





Shimadzu_FT_Right_4x4x2.mat

Load NIRS File



Generate the predicated model




SPM5

Session 1: trial specification in secs


number of conditions/trials	1
name for condition/trial 1 ?	Filger Tapping
vector of onsets - Filger Tapping	30 120 210 300 390
duration[s] (events = 0)	30

?

20 80 140 200 260
20




Filtering




The image shows a screenshot of the SPM5 software interface, specifically the 'Input for the smoothing and detrending process' dialog box. The window has a blue title bar with the SPM5 logo and standard Windows window controls (minimize, maximize, close). The main area is divided into three rows of input fields. The first row is for the 'Low-pass filter?' with a dropdown menu set to 'hrf'. The second row is for 'Detrending?' with a dropdown menu set to 'Wavelet-MDL'. The third row is for the '# of wavelet coefficients' with a text input field containing the number '4'. The background of the window is a light gray.


Input for the smoothing and detrending process	
Low-pass filter?	hrf
Detrending?	Wavelet-MDL
# of wavelet coefficients	4




Specify the Baseline



TimeSeries Analysis



Save




For example....


Shimadzu_Filtered_FT_Right_4x4x2.mat



Specify 1st Level


Specify 1st Level



NIRS Specification








SPM5

Session 1: trial specification in secs

number of conditions/trials	1
name for condition/trial 1 ?	finger tapping
vector of onsets - finger tapping	20 80 140 200 260
duration[s] (events = 0)	20

?

20 80 140 200 260




SPM5

Session 1: Other regressors

Detrending?	Wavelet-MDL
The number of wavelet coefficients at	4
Low-pass filter?	hrf
Correct for serial correlations?	<div><div>none</div><div>AR(1)</div></div>


?



Created SPM_indiv_HbO.mat

Estimate

Estimate



The image shows a Windows-style dialog box titled "NIRS_Estimation". It contains two main sections: "Estimation setup" and "Group Analysis". In the "Estimation setup" section, the "Individual Analysis" checkbox is checked, and a text field contains the file path "#shimadzu_sample#Result#SPM_indiv_HbO.mat". Below this is a button labeled "Specify SPM.mat file". The "Group Analysis" section is currently disabled, showing fields for "The number of subjects" and "The minimum number of overlapped individual subjects", along with a button "Specify the individual results". At the bottom, there are three checkboxes: "Oxy- Hb" (checked), "Deoxy- Hb" (unchecked), and "Total- Hb" (unchecked). A red oval highlights the "Estimation" button at the very bottom of the dialog.

NIRS_Estimation

Estimation setup

☒ Individual Analysis ☐ Group Analysis

#shimadzu_sample#Result#SPM_indiv_HbO.mat

Specify SPM.mat file

Group Analysis


The number of subjects

The minimum number of overlapped individual subjects

Specify the individual results

☒ Oxy- Hb ☐ Deoxy- Hb ☐ Total- Hb


Estimation





Updated SPM_indiv_HbO.mat

Results NIRS

Results NIRS









SPM contrast manager

define contrast...

name: Right Finger Tapping

type: ☒ t-contrast ☐ F-contrast

contrast weights vector: 1 0 0 0

...submit


Design matrix


parameter estimability



Reset Cancel OK


name defined, contrast not defined

?


The image shows a software window titled "SPM contrast manager". Inside, there's a section titled "define contrast...". It has three main input areas: "name" with a text box containing "Right Finger Tapping", "type" with two radio buttons ("t-contrast" is selected, "F-contrast" is not), and "contrast" with a text box containing the vector "1 0 0 0". Below the vector box is a button labeled "...submit", which is circled in red. To the right of these inputs is a large rectangular area labeled "Design matrix" showing a grayscale pattern of horizontal and vertical lines. Below the design matrix is a row of four empty boxes labeled "parameter estimability". At the bottom left of the dialog are three buttons: "Reset", "Cancel", and "OK", with the "OK" button circled in red. At the very bottom, there's a status bar that says "name defined, contrast not defined" and a small green question mark icon.









Created cinterp_SPM_nirs_HbO1_dorsal.mat



Created cinterp_SPM_nirs_HbO1_right.mat




Created cinterp_SPM_nirs_HbO1_left.mat




Created cinterp_SPM_nirs_HbO1_frontal.mat

Deoxy-Hb




T-map of **right** finger tapping



OxyHb

High t value was indicated in primary motor cortex of left side.

The result showed that increased OxyHb and decreased DeoxyHb was found at tapping




DeoxyHb

The t value of the minus was indicated in primary motor cortex of left side.

Appendix

Spatial Registration

Spatial Registration (with MRI)



Preparation

- MRI image
 - T1W
 - Normalized T1W
 - Gray Matter and White Matter images
- File for Coefficient
 - Mat file for conversion
- Location data of 4 Reference points and Probes
 - Measured using Fastrak
- Registration File
- Ch. Config File

Type	Folder name	File name
T1W	WithMRI	mri.nii
Normalized T1W	WithMRI	wmri.nii
Gray Matter	WithMRI	c1mri.nii
White Matter	WithMRI	c2mri.nii
Mat file for conversion	WithMRI	mri_sn.mat
Location data	Registration	Shimadzu_Real_4x4x2_FingerTapping.txt
Ch. Config file	Ch_Config	Shimadzu_4x4x2_48ch.txt
Saved file	Registration	Shimadzu_Saved_48ch_with_MRI

Set Reference Points and Optode

Indicator Locations

Registration of NIRS Channel With MRI

☒ 1-set of optode holder ☐ 2-set of optode holder

Real Coordinates (3D digitizer) MNI Coordinates (mm)

Reference Positions Total No.


Optode Positions Total No.

Select Real Coordinate File (xls or .txt)

Channel Positions
The spatial information of NIRS channels and optodes

MR Image Ch. Configuration NIRS-MRI Alignment View Ch. **Save Ch.**

Select Real Coordinate File



The window is titled "Indicator Locations" and contains a section for "Registration of NIRS Channel With MRI". It has two radio buttons: "1-set of optode holder" (selected) and "2-set of optode holder".

Real Coordinates (3D digitizer)

Reference Positions: Total No. 4

12.23	-0.22	10.21
25	-0.32	23.68
11.82	7.11	20.59
11.48	-8.04	20.54

Optode Positions: Total No. 32

17.19	-8.02	26.23
19.31	-9.37	23.47
19.31	-9.47	20.14
17.93	-9.22	16.99
22.26	-7.83	23.95
20.37	-7.03	26.76
21.01	-7.92	17.38
22.23	-7.95	20.59
22.98	-4.89	27.53
24.33	-5.16	23.97
24.48	-5.47	20.65
23.26	-5.53	17.57
25.38	-1.67	23.86
24.19	-1.61	27.37
24.6	-2.48	17.26
25.49	2.1	26.5
25.38	1.11	26.76

MNI Coordinates (mm)

Reference Positions


Channel Positions

The spatial information of NIRS channels and optodes


Select Real Coordinate File (xls or .txt)

MR Image Ch. Configuration NIRS-MRI Alignment View Ch. Save Ch.

MR Image (Select T1 and Normalized image)



Set Reference Points



NZ
+
CZ
+
AL
+
AR

Indicator Locations

Registration of NIRS Channel With MRI

☒ 1-set of optode holder

☐ 2-set of optode holder

Real Coordinates (3D digitizer)

Reference Positions

Total No.4

12.23	-0.22	10.21
25	-0.32	23.68
11.82	7.11	20.59
11.48	-8.04	20.54

Optode Positions

Total No.32

17.19	-8.02	26.23
19.31	-9.37	23.47
19.31	-9.47	20.14
17.93	-9.22	16.99
22.26	-7.83	23.95
20.37	-7.03	26.76
21.01	-7.92	17.38
22.23	-7.95	20.59
22.98	-4.89	27.53
24.33	-5.16	23.97
24.48	-5.47	20.65
23.26	-5.53	17.57
25.38	-1.67	23.86
24.19	-1.61	27.37
24.6	-2.48	17.26
25.49	-2.1	20.5
24.88	1.11	26.78

Select Real Coordinate File (xls or .txt)

MNI Coordinates (mm)

4 Reference Positions

0.5	90.4	-17.8
4.7	-20.9	103.3
-77	-22.1	-17.8
85.1	-12.5	-17.8

Channel Positions

The spatial information of NIRS channels and optodes

MR Image

Ch. Configuration

NIRS-MRI Alignment

View Ch.

Save Ch.

Set Ch.Config File

Indicator Locations

Registration of NIRS Channel With MRI

☐ 1-set of optode holder ☒ 2-set of optode holder

Real Coordinates (3D digitizer)

Reference Positions Total No. 4

12.23	-0.22	10.21
25	-0.32	23.68
11.82	7.11	20.59
11.48	-8.04	20.54

Optode Positions Total No. 32

17.19	-8.02	26.23
19.31	-9.37	23.47
19.31	-9.47	20.14
17.93	-9.22	16.99
22.26	-7.83	23.95
20.37	-7.03	26.76
21.01	-7.92	17.38
22.23	-7.95	20.59
22.98	-4.89	27.53
24.33	-5.16	23.97
24.48	-5.47	20.65
23.26	-5.53	17.57
25.38	-1.67	23.86
24.19	-1.61	27.37
24.6	-2.48	17.26
25.49	-2.1	20.5

Select Real Coordinate File (xls or .txt)

MNI Coordinates (mm)

4 Reference Positions

0.5	90.4	-17.8
4.7	-20.9	103.3
-77	-22.1	-17.8
85.1	-12.5	-17.8


48 Channel Positions

The spatial information of NIRS channels and optodes

MR Image **Ch. Configuration** NIRS-MRI Alignment View Ch. Save Ch.

NIRS-MRI Alignment

Select parameter file and then Gray matter and White matter file



MNI Coordinates

Indicator Locations

Registration of NIRS Channel With MRI

☐ 1-set of optode holder ☒ 2-set of optode holder

Real Coordinates (3D digitizer)

Reference Positions Total No. 4

12.23	-0.22	10.21
25	-0.32	23.68
11.82	7.11	20.59
11.48	-8.04	20.54

Optode Positions Total No. 32

17.19	-8.02	26.23
19.31	-9.37	23.47
19.31	-9.47	20.14
17.93	-9.22	16.99
22.26	-7.83	23.95
20.37	-7.03	26.76
21.01	-7.92	17.38
22.23	-7.95	20.59
22.98	-4.89	27.53
24.33	-5.16	23.97
24.48	-5.47	20.65
23.26	-5.53	17.57
25.38	-1.67	23.86
24.19	-1.61	27.37
24.6	-2.48	17.26
25.49	-2.1	20.5

Select Real Coordinate File (xls or .txt)

MNI Coordinates (mm)

4 Reference Positions

0.5	90.4	-14.2
0.5	-19.7	103.3
-78.2	-13.7	-35.8
83.9	-4.1	-35.8


48 Channel Positions

MNI Coordinates


CH01:	60	-57	37
CH02:	64	-25	32
CH03:	63	10	30
CH04:	47	-71	43
CH05:	62	-44	55
CH06:	61	-4	50
CH07:	57	24	36
CH08:	51	-56	60
CH09:	58	-24	62
CH10:	53	6	57
CH11:	40	-68	60
CH12:	51	-39	69
CH13:	52	-7	62
CH14:	49	23	52
CH15:	41	-52	71
CH16:	46	-23	72
CH17:	51	9	58
CH18:	26	-68	68

MR Image Ch. Configuration NIRS-MRI Alignment **View Ch.** Save Ch.

View Ch.




Save Ch.



Shimadzu_Saved_48ch_with_MRI.mat

Spatial Registration(without MRI)



Preperation

- Location data of Reference point
 - Shimadzu_4x4x2_FingerTapping_origin.csv
- Location data of Probe and Channel
 - Shimadzu_4x4x2_FingerTapping_others.csv

rain.csv

Microsoft Excel - 4x4x2_FingerTapping_origin.csv

ファイル(F) 編集(E) 表示(V) 挿入(I) 書式(O) ツール(T) データ(D) ウィンドウ(W) ヘルプ(H)

MS Pゴシック 11 B I U

校閲結果の返信(C)... 校閲結果の差し込み終了(N)...

	A1		fx	Label					
	A	B	C	D	E	F	G	H	I
1	Label	X	Y	Z					
2	NzHS	12.23	-0.22	10.21					
3	IzHS								
4	ARHS	11.48	-8.04	20.54					
5	ALHS	11.82	7.11	20.59					
6	Fp1 HS								
7	Fp2 HS								
8	FzHS								
9	F3HS								
10	F4HS								
11	F7HS								
12	F8HS								
13	CzHS	25	-0.32	23.68					
14	C3HS								
15	C4HS								
16	T3HS								
17	T4HS								
18	PzHS								
19	P3HS								
20	P4HS								
21	T5HS								
22	T6HS								
23	O1 HS								
24	O2HS								
25									
26									

4x4x2_FingerTapping_origin/

図形の調整(R) オートシェイプ(U)

コマンド

Shimadzu_4x4x2_FingerTapping_others.csv

	A	B	C	D	E	F	G	H	I
14	R7	24.19	-1.61	27.37					
15	T8	24.6	-2.48	17.26					
16	R8	25.49	-2.1	20.5					
17	T9	24.28	1.11	26.78					
18	R9	25.61	0.68	23.3					
19	T10	25.52	0.59	19.92					
20	R10	24.78	0.27	16.65					
21	T11	24.84	3.89	22.73					
22	R11	23.67	4.52	26.55					
23	T12	23.96	3.04	16.17					
24	R12	24.94	4	19.52					
25	T13	21.18	6.6	25.62					
26	R13	22.74	6.14	22.31					
27	T14	22.46	5.91	18.83					
28	R14	21.36	5.41	15.77					
29	T15	19.99	8.09	21.72					
30	R15	18.35	7.27	24.85					
31	T16	18.38	6.79	15.5					
32	R16	19.7	7.34	18.52					
33	CH01	18.25	-8.695	24.85					
34	CH02	19.31	-9.42	21.805					
35	CH03	18.62	-9.345	18.565					
36	CH04	18.78	-7.525	26.495					
37	CH05	20.785	-8.6	23.71					
38	CH06	20.77	-8.71	20.365					
39	CH07	19.47	-8.57	17.185					

T**, X, Y, Z

R**, X, Y, Z

CH**, X, Y, Z

CH is center
point of logical
target T, R

NIRS_Registration_Standalone

NIRS_Registration_Standalone

Registration of NIRS Channel Without MRI

☐ Without 3D Digitizer ☒ With 3D Digitizer

☐ 1-set of optode holder ☒ 2-set of optode holder

Input panel of MNI coordinates

Input point is about ☒ Optodes ☐ Ch. Config. ☐ Channels

Select the file to contain MNI coordinates of NIRS optodes

Select the SPM template to specify MNI coordinates

MNI coordinate of specific optodes:[Optd #, x, y, z]

Add Delete Clear all

From Real Coordinates to MNI Space

Select the file (Reference position in REAL space)

Select the file (Optode and Channel position in REAL space)

Registration (use the NFRI function)

Project MNI Coordinate to Rendered Brain

View the Optodes on the rendered brain

保存

The spatial information of NIRS channels and optodes

If you have the 3D digitizing system,
NFRI (NIRS tools (Singh et al. 2005) incorporated in NIRS-SPM
allow spatial registration of NIRS channels from real space to
MNI space.
In order to produce results using this function, you are also
required to cite the following paper in addition to NIRS-SPM
papers (Ye et al., 2009; Jang et al., 2009):
A.K. Singh et al. / NeuroImage 27 (2005) 842-851

NIRS_Registration_Standalone

Registration of NIRS Channel Without MRI

☐ Without 3D Digitizer ☒ With 3D Digitizer
☐ 1-set of optode holder ☒ 2-set of optode holder

Input panel of MNI coordinates

Input point is about ☒ Optodes Ch. Config. ☐ Channels

Select the file to contain MNI coordinates of NIRS optodes

Select the SPM template to specify MNI coordinates

MNI coordinate of specific optodes: [Optd #, x, y, z]

 Add Delete Clear all

From Real Coordinates to MNI Space

Select the file (Reference position in REAL space)

4x4x2_FingerTapping_origin.csv

Select the file (Optode and Channel position in REAL space)

4x4x2_FingerTapping_others.csv

Registration (use the NFRI function)

Project MNI Coordinate to Rendered Brain


View the Optodes on the rendered brain


保存

The spatial information of NIRS channels and optodes


32 Optode Position and 48 Channel Position

MNI Estimation







Estimation result




Estimation result




Zoomed



Project....



Save Ch.



Shimadzu_Saved_48ch_without_MRI.mat